

Intercomparison of Mobility Particle Size Spectrometers

Project No.: MPSS-2017-5-4

Principal Investigator: Noemi Perez

Home Institution: IDAEA-CSIC
Jordi Girona 18-26
08034 Barcelona
Spain

Participant: Noemí Perez

Candidate: **ES-IDAEA-CSIC**
Made by: **TROPOS Home-Made**
Counter (SN): TSI CPC Model 3772, SN: 3772140203
Software: TROPOS

Location of the quality assurance: TROPOS Leipzig, lab 118

Comparison period: October 09, 2017 – October 13, 2017

Last Intercomparison (with Project No.):

Summary of Intercomparison:*Pre-Status:*

The instrument arrived with participant. During the Pre-Status, the performance of the system showed a concentration 11% lower than the TROPOS Reference Instrument No.4. The PSL check showed a correct peak at 204.99 nm. The TROPOS software version should be updated to the latest version 6.66; therefore, the whole PC with Windows XP has to be changed. In 2018, ES-IDAEA-CSIC will get an offer from TROPOS to update the hardware and a possibility to change the instrument to an electronic box with USB-NI connection. The pressure drop sensor in the aerosol inlet is shifting.


Final-Status:

During the Final-Status, the performance of the system showed a concentration 1% higher than the TROPOS Reference Instrument No.4. The candidate was operating with a cleaned DMA and settings were updated in the TROPOS evaluation program. During the whole workshop, ES-IDAEA-CSIC used a Kr.85 source from TROPOS. The candidate passed the quality standards of ACTRIS and GAW.

Information about the instruments:**Date of check: October 09, 2017**

<i>List of Components</i>	TROPOS Reference MPSS No.1	TROPOS Reference MPSS No.4	Candidate
<i>Position</i>	Line 1	Line 2	Line 2
<i>Company</i>	TROPOS	TROPOS	TROPOS
<i>Software</i>	TROPOS	TROPOS	TROPOS V4.7.2 (win xp)
<i>CPC-MPSS</i>	TSI CPC, Model 3772	TSI CPC, Model 3772	TSI CPC, Model 3772
<i>CPC-total</i>	TSI CPC, Model 3010	TSI CPC, Model 3010	-
<i>flow ratio</i>	1.0 : 5.0	1.0 : 5.0	1.0 : 5.0
<i>source</i>	Kr85	Ni63	Kr85 from TROPOS
<i>HV power supply</i>	Positive	Positive	Positive
<i>DMA</i>	Hauke medium	Hauke medium	Hauke medium
<i>aerosol dryer</i>	✓	✓	✓
<i>aerosol RH- sensor</i>	✓	✓	✓
<i>aerosol T-sensor</i>	✓	✓	✓
<i>sheath RH-sensor</i>	✓	✓	✓
<i>sheath T-sensor</i>	✓	✓	✓
<i>Sheath dryer</i>	✓	✓	✓
<i>pressure sensor</i>	✓	✓	-
<i>info</i>			Offer for update MPSS

Laboratory setup:



Status of the instruments:

Date of check (Pre-Status): October 09, 2017

<i>CPC status</i>	MPSS		Total CPC	
<i>power/status</i>	LED green	-	-	-
<i>saturator temp</i>	39.0	°C	-	°C
<i>condenser temp</i>	22.0	°C	-	°C
<i>optics temp</i>	40.0	°C	-	°C
<i>cabinet temp</i>	32.4	°C	-	°C
<i>ambient pressure</i>	99.3	kPa	-	kPa
<i>orifice pressure</i>	80.8	kPa	-	kPa
<i>nozzle pressure</i>	2.6	kPa	-	kPa
<i>laser current</i>	39	mA	-	mA
<i>liquid level</i>	full	-	-	-

Date of check (Final-Status): October 12, 2017

<i>CPC status</i>	MPSS		Total CPC	
<i>power/status</i>	LED green	-	-	-
<i>saturator temp</i>	39.0	°C	-	°C
<i>condenser temp</i>	22.0	°C	-	°C
<i>optics temp</i>	40.0	°C	-	°C
<i>cabinet temp</i>	34.4	°C	-	°C
<i>ambient pressure</i>	99.7	kPa	-	kPa
<i>orifice pressure</i>	82.7	kPa	-	kPa
<i>nozzle pressure</i>	2.7	kPa	-	kPa
<i>laser current</i>	40	mA	-	mA
<i>liquid level</i>	full	-	-	-

Date of system checks:

<i>date</i>	09.10.2017	10.10.2017	11.10.2017	13.10.2017	unit
<i>total CPC flow</i>	-	-	-	-	l/min
<i>aerosol flow (DMA)</i>	-	-	-	-	l/min
<i>aerosol flow (UDMA)</i>	-	-	-	-	l/min
<i>aerosol flow (total)</i>	0.985	-	1.006	1.008	l/min
<i>Zero MPSS</i>	0	-	0	0	#/cm ³
<i>Zero total CPC</i>	-	-	-	-	#/cm ³
<i>PSL 203 nm</i>	204	-	204	201	nm
<i>HV – 0 V</i>	0.2	-	0.3	-	V
<i>HV – 5 V</i>	5.1	-	5.4	-	V
<i>HV – 100 V</i>	99.7	-	99.8	-	V
<i>HV – 1000 V</i>	999.8	-	1000.7	-	V

Special Information regarding the Candidate:

<i>Was it necessary to:</i>	yes/no (date)	old part (ID/SN)	new part (ID/SN)	information
<i>clean the aerosol inlet</i>	Yes	-	-	replace capillary
<i>change aerosol Nafion dryer</i>	Yes	MT052913-05-3	MT112916-01-12	to old
<i>change sheath Nafion dryer</i>	Yes	ND0.7-05a	ND0.7-26e	to old
<i>check source</i>	Yes	-	-	TROPOS source
<i>change HV power supply</i>	No	-	-	-
<i>clean/change DMA</i>	Yes	-	-	Cleaned; DMA okay
<i>change aerosol RH/T-sensor</i>	No	-	-	-
<i>change sheath RH/T-sensor</i>	No	-	-	-
<i>change pressure sensor</i>	No	-	-	-
<i>change inlet Nafion dryer (500)</i>	No	-	-	-
<i>Change Total filter</i>	No	-	-	-

PSL Scan and calibration: Latex 203 nm +/- 4 nm




Figure 01: Measurement of latex 203 nm: Particle size distribution (raw concentration) for latex 203 nm on October 11th, 2017.




Figure 02: Measurement of latex 203 nm: Particle size distribution (raw concentration) for latex 203 nm on October 11th, 2017.

Status of the TROPOS Reference Instruments: Particle Number Size Distribution




Figure 03: Comparison of mean particle number size distribution of TROPOS Reference MPSS No.1 against TROPOS Reference MPSS No.4 from October 09, 2017 08:00 PM – October 10, 2017 06:00 AM. Multiple charge correction, internal diffusion losses and CPC efficiency are included.

Status of the TROPOS Reference Instruments: Time Series




Figure 04: Time series (October 09, 2017 08:00 PM – October 10, 2017 06:00 AM) of the integrated particle number concentration ($N_{10-800nm}$) of the MPSS and total number concentration (N_{total}) of the Reference TSI-CPC Model 3010. Multiple charge correction, internal diffusion losses and CPC flow corrections are included.

Status of the TROPOS Reference Instruments Correlation




Figure 05: Linear regression between the number concentrations of the TROPOS total Reference TSI CPC Model 3010 No.1 and TROPOS total Reference TSI CPC Model 3010 No.4. Coincidence corrections and CPC flow corrections are included.




Figure 06: Linear regression between the number concentrations of the TROPOS Reference MPSS No.1 and TROPOS Reference MPSS No.4. Multiple charge correction, internal diffusion losses and CPC flow corrections are included.

Pre-Status of the Candidate: Particle Number Size Distribution




Figure 07: Comparison of mean particle number size distribution of TROPOS Reference MPSS No.4 against ES-IDAEA-CSIC from October 09, 2017 08:00 PM – October 10, 2017 06:00 AM. Multiple charge correction, internal diffusion losses and CPC efficiency are included.

Pre-Status of the Candidate: Time Series




Figure 08: Time series (October 09, 2017 08:00 PM – October 10, 2017 06:00 AM) of the integrated particle number concentration ($N_{10-800nm}$) of the MPSS and total number concentration (N_{total}) of the Reference TSI-CPC Model 3010. Multiple charge correction, internal diffusion losses and CPC flow corrections are included.

Pre-Status of the Candidate: Correlation




Figure 09: Linear regression between the number concentrations of the TROPOS Reference MPSS No.4 and ES-IDAEA-CSIC. Multiple charge correction, internal diffusion losses and CPC flow corrections are included.

Final-Status of the Candidate: Particle Number Size Distribution




Figure 10: Comparison of mean particle number size distribution of TROPOS Reference MPSS No.4 against ES-IDAEA-CSIC from October 13, 2017 06:00 PM – October 16, 2017 06:00 AM. Multiple charge correction, internal diffusion losses and CPC efficiency are included.

Final-Status of the Candidate: Time Series




Figure 11: Time series (October 13, 2017 06:00 PM – October 16, 2017 06:00 AM) of the integrated particle number concentration ($N_{10-800\text{nm}}$) of the MPSS and total number concentration (N_{total}) of the Reference TSI-CPC Model 3010. Multiple charge correction, internal diffusion losses and CPC flow corrections are included.

Final-Status of the Candidate: Correlation




Figure 12: Linear regression between the number concentrations of the TROPOS Reference TSI CPC Model 3010 SN: 2337 and TROPOS Reference MPSS No.4 (October 13, 2017 06:00 PM – October 16, 2017 06:00 AM). Multiple charge correction, internal diffusion losses and CPC flow corrections are included.




Figure 13: Linear regression between the number concentrations of the TROPOS Reference TSI CPC Model 3010 SN: 2337 and ES-IDA EA-CSIC (October 13, 2017 06:00 PM – October 16, 2017 06:00 AM). Multiple charge correction, internal diffusion losses and CPC flow corrections are included.




Figure 14: Linear regression between the number concentrations of the TROPOS Reference MPSS No.4 and ES-IDA EA-CSIC (October 13, 2017 06:00 PM – October 16, 2017 06:00 AM). Multiple charge correction, internal diffusion losses and CPC flow corrections are included.