

Leibniz-Institut für Troposphärenforschung Permoserstraße 15 04318 Leipzig



Leibniz Institute for Tropospheric Research

CPC Model: TSI CPC 3750

CPC Serial Number: 3750193902

Customer: TSI Instruments Ltd.

Description: Calibration of a Condensation Particle Counter (CPC, Model 3750)

Date of Calibration: November 13, 2019

Summary of Intercomparison:

The candidate passed the quality standards of ACTRIS and GAW. The candidate reached 99% efficiency at 40 nm. The Dp50 is at 6.57 nm. The CPC efficiency curve corresponds to the standard of ACTRIS and GAW.

Certificate / Reference: WCCAP


Date of issue: November 13, 2019 Signature:

Reviewed by: TROPOS Name: Kay Weinhold

Page 1 / 4

Leibniz-Gemeinschaft







Leibniz Institute for Tropospheric Research

Date of arrival of instrument in calibration lab: November 11, 2019

Instrument: Condensation Particle Counter Model and serial number of instrument: CPC 3750 S/N 3750193902

Result of physical inspection: no damages

Result of functional test: functional test successful, no problems

nominal flow rate 1.0 l/min Internal parameters of instrument

Model and identification number of

aerosol electrometer: TSI Electrometer Model 3068, S/N 70838596

Electrometer calibration certificate: September 5, 2018, calibrated at PTB

Braunschweig

Corrections of electrometer, for instance,

differing flow rate:

Within tolerance range (+/-2%); reference: 4.0

I/min, measured: 4.000 I/min

Software for recording: LabView 2010; National Instruments; Program

"LabCount.vi"

Date of calibration: November 13, 2019 Lab temperature and pressure: 24.5°C, 982.5 mbar

Measured aerosol flow rate of CPC: 1.007 I/min

Uncertainty in measured flow rate: 3%

Flowmeter used: Gilian Gilibrator V; S/N 1711008-S,

January, 2018

Particles and gases used for calibration: silver particles and nitrogen Method of particle generation: tube furnace generator Zero measurement of instrument: 0 particles/cm3 in 5 minutes


Results (using pulse output):

Particle size (nm)	40	30	20	15	10		
Number concentration (cm-3)	1277	1270	1512	1275	1290		
Counting efficiency η	0.99	0.99	1.00	1.00	0.90		
Particle size (nm)	09	08	07	06			
Number concentration (cm-3)	1371	876	814	510			
Counting efficiency η	0.82	0.72	0.57	0.37			

SWIFT CODE: COBADEFF 860









Leibniz Institute for Tropospheric Research




Fig. 1: Counting efficiency for CPC 3750 S/N 3750193902 against aerosol electrometer 3068 S/N 70838596; silver particles between 6 and 40 nm were used for calibration; the calculated Dp50 is 6.57 nm.

Status information:

Status	T SAT	T CON	T OPT	T CAB	P AMB	P VAC
from display	39.0	18.0	40.0	24.3	99.0	79.2
Status	P OR	P NO	Laser	LV	flow	P INLET
from display	76.4	2.33	39	full	1.007	-0.4

Results (using pulse output):

Results (using pulse output):							
without coincidence correction							
Concentration EM in #/cm³	66849	58164	50013	43197	32294		
Number concentration without							
coincidence correction (cm-3)	52691	46903	41157	36249	28000		
Counting efficiency η	0.79	0.80	0.82	0.84	0.86		
Concentration EM in #/cm³	19887	11244	5127	1177			
Number concentration without							
coincidence correction (cm-3)	17978	10519	4956	1159			
Counting efficiency η	0.90	0.93	0.96	0.98			

Page 3 / 4



Commerzbank Leipzig KTO 102 14 50 BLZ 860 400 00 IBAN: DE77 8604 0000 0102 1450 00 SWIFT CODE: COBADEFF 860 Mitglied der Leibniz-Gemeinschaft







Leibniz Institute for Tropospheric Research

with coincidence correction							
Concentration EM in #/cm³	66849	58164	50013	43197	32294		
Number concentration with	60566	E0222	F0627	105.10	22226		
coincidence correction (cm-3)	68566	59332	50627	43542	32326		
Counting efficiency η	1.02	1.02	1.01	1.00	1.00		
Concentration EM in #/cm³	19887	11244	5127	1177			
Number concentration with							
coincidence correction (cm-3)	19907	11349	5265	1235			
Counting efficiency η	1.00	1.00	1.02	1.04			




Fig. 2: Linearity test for TSI CPC 3750 SN 3750193902 against aerosol electrometer 3068 SN 70838596; silver particles with a diameter of 30 nm were used for number concentrations between 1000 and 70000 particles per cm³.

Date of issue: November 13, 2019

Reference: TSI electrometer, model 3068, SN 70838596

Reviewed: TROPOS / Kay Weinhold

Page 4 / 4