

Leibniz-Institut für Troposphärenforschung Permoserstraße 15 04318 Leipzig



Leibniz Institute for Tropospheric Research

CPC Model:	TSI CPC 3772
------------	--------------

CPC Serial Number: 71011008

Customer: -

**Description:** Calibration of a Condensation Particle Counter (CPC, Model 3772)

**Date of Calibration:** June 04, 2020

#### **Summary of Intercomparison:**

The candidate passed the quality standards of ACTRIS and GAW. The candidate reached 100% efficiency at 40 nm. The Dp50 is at 9.89 nm. The candidate was calibrated to Dp50 of 10nm. The CPC efficiency curve corresponds to the standard of ACTRIS and GAW.

Certificate / Reference: WCCAP

Date of issue: June 04, 2020 Signature:

Reviewed by: TROPOS Name: Kay Weinhold

Page 1 / 4

Leibniz-Gemeinschaft







### Leibniz Institute for Tropospheric Research

Date of arrival of instrument in calibration lab: February 11, 2020

Instrument: Condensation Particle Counter Model and serial number of instrument: CPC 3772 S/N 71011008

Result of physical inspection: no damages

Result of functional test: functional test successful, no problems

Internal parameters of instrument nominal flow rate 1.0 l/min

Model and identification number of

aerosol electrometer: TSI Electrometer Model 3068, S/N 70838596

**Electrometer calibration certificate:** September 5, 2018, calibrated at PTB

Braunschweig

Corrections of electrometer, for instance,

differing flow rate: Within tolerance range (+/-2%); reference: 4.0

I/min, measured: 4.000 I/min

Software for recording: LabView 2010; National Instruments; Program

"LabCount.vi"

Date of calibration: June 04, 2020 Lab temperature and pressure: 23.0°C, 982.0 mbar

Measured aerosol flow rate of CPC: 1.020 I/min

Uncertainty in measured flow rate: 3%

Flowmeter used: Gilian Gilibrator V; S/N 1711008-S,

January, 2018

Particles and gases used for calibration: silver particles and nitrogen Method of particle generation: tube furnace generator Zero measurement of instrument: 0 particles/cm3 in 5 minutes

Results (using pulse output): Pre-Status


Particle size (nm)	40	30	20	15	10
Number concentration (cm-3)	1029	1374	1182	1451	1237
Counting efficiency η	1.01	1.00	1.01	0.99	0.85
Particle size (nm)	09				
Number concentration (cm-3)	1079				
Counting efficiency η	0.77				

Page 2 / 4



SWIFT CODE: COBADEFF 860







## Leibniz Institute for Tropospheric Research

#### Results (using pulse output): After calibrating

Particle size (nm)	40	30	20	15	10
Number concentration (cm-3)	1429	1243	1679	1319	775
Counting efficiency η	1.01	1.01	0.97	0.88	0.52
Particle size (nm)	09	08	07		
Number concentration (cm-3)	567	212	15		
Counting efficiency η	0.37	0.15	0.01		





Fig. 1: Counting efficiency for CPC 3772 S/N 71011008 against aerosol electrometer 3068 S/N 70838596; silver particles between 7 and 40 nm were used for calibration; The instrument was calibrated to Dp50 of 10nm and resulted in a Dp50 of 9.89nm. The graph shows the counting efficiency before the calibration and after.

Commerzbank Leipzig KTO 102 14 50 BLZ 860 400 00 IBAN: DE77 8604 0000 0102 1450 00 SWIFT CODE: COBADEFF 860

Leibniz-Gemeinschaft







# Leibniz Institute for Tropospheric Research

#### **Status information:**

Status	T SAT	T CON	T OPT	T CAB	P AMB	P VAC
from display	39.0	25.5	40.0	29.8	98.1	-
Status	P OR	P NO	Laser	LV	flow	P INLET
from display	75.6	2.6	50	full	1.020	-

Date of issue: June 04, 2020

Reference: TSI electrometer, model 3068, SN 70838596

Reviewed: TROPOS / Kay Weinhold

Page 4 / 4

http://www.tropos.de

Commerzbank Leipzig KTO 102 14 50 BLZ 860 400 00 IBAN: DE77 8604 0000 0102 1450 00 SWIFT CODE: COBADEFF 860 Mitglied der Leibniz-Gemeinschaft